



Prescribed for Progress™

SXT7000e² Medical

EASY ENTRY UPRIGHT ELLIPTICAL

With all the same total body exercise functionality as the SXT7000e², this model has been certified for use as a Medical Device in hospitals and rehab clinics.

SXT7000e² Medical Certification



- Medical CE Class IIa Certification BS EN ISO13485:2008
- Calibrated for medical accuracy
- Delivers exact training protocols for watts (constant work) and RPM
- Full-length handrails and 50% lower step-up platform for easy-entry
- Bio-Flex™ footbeds improve foot and ankle circulation
- Pedal movement simulates uneven surface to recruit stabilizing muscles
- True-to-life, natural walking motion
- Optimised torso rotation and core muscle recruitment for a weight-bearing total body workout
- Safe and easy entry for progression to weight-bearing elliptical motion
- Bi-directional exercise
- Very low starting resistance

SPECIFICATIONS

Specification	SXT7000e ² Medical
Workout Position	Upright
Console	Intelli-Fit
Seat System	n/a
Hangrips & Footbeds	Dual-position handgrips and orthopedic footbeds with patented Bio-Flex™ technology
Easy Entry Step & Handrails	Yes
Medical CE IIa Certification	Yes
Contact Heart Rate	Yes
Resistance	Workload levels range from 1 - 20, adjustable in 0.1 increments, providing 191 levels of resistance
Programs	Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training
Feedback	Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS
Power	Self-generating from just 6 watts with auto recharge and battery back up
Length	190cm
Width	81cm
Height	160cm
Weight	138kg
Max. User Capacity	193kg
Color	Silver powder coated frame and cool grey covers
Warranty	3 years parts

CE 0434

Medical Device CE Class IIa certification BS EN ISO13485:2008 for use in hospitals and clinics where calibrated accuracy is required for watts and RPM.

