


## ELITE TILT TABLE

The ELITE Tilt Table is a Variable Height Tilt Table with hand switch controlled actuators. As well as the standard ELITE rigid design it offers as standard swing away, adjustable, split footplates and split leg section for individual securing of the legs. Individual braked castors provide ease of movement as well as a stable base throughout the treatment session. The base design allows for easy access with patient lifting devices and the variable height allows for easy transfer from wheelchair or plinth. A full tilt range from 0 through 90 degrees is standard.

Includes activity tray (arm rests are optional).



### ***SPECIFICATIONS:***

Minimum Height	50cm
Maximum Height	85cm
Cushion Width	65cm
Cushion Length	190cm
Foot Plate width	45cm
Centre Section	67cm
Minimum Tilt	0° (horizontal)
Maximum Tilt	90° (vertical)
Total Length	196cm
Retractable Wheels	Standard

### VINYL COLOURS:

RAVEN  
BLACK

STORM  
GRAY

SKY  
BLUE

OCEAN  
BLUE

BALMAIN  
BLUE

WADI  
BEIGE



**metron**  
advanced health technology

57 Aster Avenue, Carrum Downs VIC 3201, Australia  
Phone: +61 (0) 3 9775 1234 Fax: +61 (0) 3 9775 1990